FOR PROFESSIONAL DENTAL USE ONLY.

Clean and sterilize the applicator prior to use (see 2. Cleaning and sterilization).

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For use with capsulated dental materials.

### 1. Step-by-step instructions

- Hold the capsule firmly against a hard, flat surface with the plunger facing down.
- Push down firmly, forcing the plunger into the body of the capsule until it stops.
- Mix the capsule in a high-frequency oscillating or rotating capsule mixer.
- Insert the capsule into the applicator using the grooves to hold it in place.
- Gently squeeze the dosing lever two times to push the material towards the tip.
- Extrude the material from the capsule by pressing the lever an additional two times.
- After application, remove the empty capsule using the applicator release button.

### 2. Cleaning and sterilization

#### 2.1. Manual cleaning

Step	Solution	Time (minutes)	Temperature	Instruction
1	Water	2	Ambient	Rinse soiled applicator under ambient tap water. Use a soft bristled brush if needed.
2	Hospital Grade Enzymatic detergent	20	Room temperature	Immerse and soak for required time. Clean thoroughly - scrub all external surfaces with a soft bristle brush until all visible soil has been removed. It is important to make sure that all joints and handles are effectively cleaned. Use a small diameter brush or pipe cleaner if needed. Inspect for visible soil on exposed surfaces.
3	Water	2	Warm, as delivered from hot tap	Rinse applicator under warm tap water for required time.
4	Air	As required	Ambient	Allow to air dry in clean area. Blow lumens with clean air using filtered air source or syringe.

#### 2.2. Mechanical cleaning

The applicator can be mechanically cleaned. However, due to several numbers of brands and quality of automated cleaners on the market, it is not possible to make a claim in compatibility with washer disinfectant meeting all requirements of ISO 15883 series.

The applicator should be positioned to ensure surfaces are exposed to the action of the water spray without dislodgement, and to allow drainage.

For automated cleaning of the applicator it is recommended to use an approved medical washer disinfectant. Follow instrument washer manufacturer's instructions for use.

The applicator should be dry before sterilization.

#### 2.3. Sterilization

A steam sterilizer is approved for the applicator. Proceed with sterilization procedure according to autoclave manufacturer's instructions for use.

Maximum sterilization temperature: 134°C (273.2 °F).

#### 2.3.1. Recommended exposure times:

Pre-Vac pouched cycle type. 30 minutes at 121°C (249.8 °F). 20 minutes dry time.

Pre-Vac pouched cycle type. 4 minutes at 132°C (269.6 °F). 20 minutes dry time.

### 3. Maintenance

Regularly apply dental instrument lubricant between applicator handles.

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#### Symbols

	Batch code	Rx ONLY
	Catalogue number	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Consult instructions for use	
	Manufacturer	
	Caution	
	Date of manufacture	

#### Limitation of Liability

If this Product is found defective within the Warranty Period, your exclusive remedy and Doxa Dental AB's sole obligation and liability shall be the repair or replacement of the Doxa Dental AB Product. EXCEPT WHERE PROHIBITED BY LAW, DOXA DENTAL AB WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGE ARISING FROM THE PRODUCT, WHETHER DIRECT, INDIRECT, SPECIAL, PUNITIVE, INCIDENTAL OR CONSEQUENTIAL INCLUDING LOSS OF PROFITS, REGARDLESS OF THE THEORY ASSERTED, INCLUDING WARRANTY, BREACH OF CONTRACT, NEGLIGENCE OR STRICT LIABILITY.

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Ceramir® er et varemærke tilhørende Doxa Dental AB.
Ceramir® on Doxa Dental AB:n tavaramerkki.